




Chibi 0.3(Comma san)
 96 WJacket Tubes 0.3ml Type with lid lock and anti-rotation function (working volum 0.3ml)



- Polypropylene
- Autoclavable
- DMSO compatible

BIFile CERTIFIED RNase DNase SAFE
 The product is certified RNase and DNase Free.



• 96 formatted WJacket Tubes 0.3ml Type for External Screw cap (working volume 0.3ml)

Product code	Specification	unit	retail price
WJRDS-0XS(EX)	without cap, bulk	960 pcs	JPY 57,750
WJRDS-2XS(EX)UL	without cap, Racked, Rack ID : Laser etching	96 pcs×10 racks	JPY 63,000
JD-9RT-XS(EX)UL	Empty Rack, Rack ID : Laser etching	10 racks	JPY 9,400
JD-9RT-XS(EX)	Empty Rack, Rack ID : N/A	10 racks	JPY 9,000

* When ordering, please enter the color code of the cap shown right below in the □ part of the product number.

• Screw cap for 96 format



- Adopts a screw cap that does not use O-rings, packing, or gaskets.
- Since no sealing material is used, it has good sealing properties and is safe to use.

Product code	Specification	unit	retail price
JSC-OS(EX) (□)	bulk cap	960 pcs	JPY 7,500

* When ordering, please enter the color code in the □ at the end of the product number. (Standard color is R. Example : [JSC-OS(EX) (R)])

Color option for screw cap




Function of Rack ID

- The same alphanumeric information is laser engraved on the bottom of the rack as a 2D code and on the side as a 1D barcode.
- Code specifications [2D code on the bottom (DataMatrix ECC200) + 1D barcode on the side (Code 128) (2-digit alphanumeric characters, 8-digit numbers)]
- By attaching 2D code for Rack ID onto the bottom of Rack, Tube ID and Rack ID can be read at the same time, reducing stress on the sample during reading process. (*However, only when using a reader that can read all at once)

Rack ID side engraving section

Rack ID bottom engraving section






*The image shows a 2.0 ml type product.

Chibi 0.3(Comma san)
96 WJacket Tubes 0.3ml Type with lid lock and anti-rotation function (working volum 0.3ml)



- Polypropylene
- Autoclavable
- DMSO compatible

CERTIFIED RNase DNase SAFE
The product is certified RNase and DNase Free.



• 96 formatted WJacket Tubes 0.3ml Type for External Screw cap (working volume 0.3ml)

Product code	Specification	unit	retail price
WJRDS(□)-5XS(EX)UL	Capped, Racked, Rack ID : Laser etching, Non-sterile	96 pcs × 10 racks	JPY 75,000
WJRDS(□)-6XS(EX)	Capped, bulk Non-sterile	96 pcs × 10 bags	JPY 69,000
WJRDS(□)-5XS(EX)ULS	Capped, Racked, Rack ID : Laser etching, Sterile	96 pcs × 10 racks	JPY 81,000
WJRDS(□)-6XS(EX)S	Capped, bulk Sterile	96 pcs × 10 bags	JPY 75,000

* When ordering, please enter the color code of the cap shown right below in the □ part of the product number.

• Screw cap for 96 format



- Adopts a screw cap that does not use O-rings, packing, or gaskets.
- Since no sealing material is used, it has good sealing properties and is safe to use.

Product code	Specification	unit	retail price
JSC-OS(EX)(□)	bulk cap	960 pcs	JPY 7,500

* When ordering, please enter the color code in the □ at the end of the product number. (Standard color is R. Example : [JSC-OS(EX)(R)])




Function of Rack ID

- The same alphanumeric information is laser engraved on the bottom of the rack as a 2D code and on the side as a 1D barcode.
- Code specifications [2D code on the bottom (DataMatrix ECC200) + 1D barcode on the side (Code 128) (2-digit alphanumeric characters, 8-digit numbers)]
- By attaching 2D code for Rack ID onto the bottom of Rack, Tube ID and Rack ID can be read at the same time, reducing stress on the sample during reading process. (*However, only when using a reader that can read all at once!)

Rack ID side engraving section

Rack ID bottom engraving section



*The image shows a 2.0 ml type product.